

Certificate of Analysis

Product Fetal Bovine Serum Collected in South America, sterile-filtered

Cat. No. AB-FBS-0500, AB-FBS-1050

Lot. No. AB20214240 Expiry Date 01/2026
Origin South America Manufacturing Date 01/2021

Storage ≤-15°C

Quality Profile	Specification	Result	Unit
Appearance	Clear amber	Pass	
рН	6.8 – 8.2	6.90	
Osmolality ⁽¹⁾	280 – 340	334	mOsm/kg
Endotoxin ⁽²⁾	As reported	3.0	EU/ml
Hemoglobin ⁽³⁾	≤ 25	24.10	mg/dl
Glucose ⁽⁴⁾	As reported	4.10	mM
Calcium ⁽³⁾	As reported	3.28	mM
Phosphorous ⁽³⁾	As reported	3.59	mM
GOT/AST ⁽³⁾	As reported	41	U/I
LDH ⁽³⁾	As reported	514	U/I

Protein Profile			
Total Protein ⁽⁵⁾	3.0 – 4.5	3.6	g/dl
Albumin ⁽⁶⁾	As reported	1.8	g/dl
IgG ⁽⁷⁾	≤ 250	80.12	mg/l

Sterility Testing		
Bacteria and Fungi®	Not detected	Not detected
Mycoplasma ⁽⁹⁾	Not detected	Not detected

Virus Testing ⁽¹⁰⁾		
Bovine Herpes Virus (BHV-1 and BHV-5)	Not detected	Not detected
Pestivirus (e.g. BVD-V)	Not detected	Not detected
Parainfluenza Virus Type 3 (PI-3)	Not detected	Not detected

Biological Performance			
Cell Growth MDCK (Canine Kidney Epithelial Cell Line)(11)	≥ 75	99	%
Cell Growth CHO (Chinese Hamster Ovary Cell Line)(11)	≥ 75	98	%

References/Methods

(3)

) Freeze Point Depression (7) ELIS

(2) LAL-Test (8) Membrane Filtration (total plate count, aerobic and anaerobic bacteria, yeasts and fungi)

Photometry (9) PC

(4) Hexokinase Assay (10) Cytopathic Effect in Cell Culture and Immunofluorescence

(5) Biuret Colorimetry(6) Nephelometry(11) Cell Promoting Assay

Intended Use For research use only.

Date June 25, 2021

This document has been created digitally and is valid without a signature.