

Certificate of Analysis (Preliminary)

Product FBS Xtra, Fetal Bovine Serum, Collected in South America, sterile-filtered

Cat. No. AB-FBS-0500, AB-FBS-1050

Lot. No. **Expiry Date** 06/2024 AB20214217 Origin Manufacturing Date 06/2021 South America

Storage ≤-15°C

Quality Profile	Specification	Result	Unit
Appearance	Clear amber	Pass	
рН	6.8 – 8.2	7.28	
Osmolality ⁽¹⁾	280 – 340	321	mOsm/kg
Endotoxin ⁽²⁾	As reported	< 0.36	EU/mI
Hemoglobin ⁽³⁾	≤ 25	9.3	mg/dl
Glucose ⁽⁴⁾	As reported	1.11	mM
Calcium ⁽³⁾	As reported	1.65	mM
Phosphorous ⁽³⁾	As reported	5.94	mM
GOT/AST ⁽³⁾	As reported	29	U/I
LDH ⁽³⁾	As reported	326	U/I

Protein Profile			
Total Protein ⁽⁵⁾	3.0 – 4.5	2.2	g/dl
Albumin ⁽⁶⁾	As reported	1.3	g/dl

Sterility Testing			
Bacteria and Fungi®	Not detected	Not detected	
Mycoplasma ⁽⁹⁾	Not detected	Not detected	

Virus Testing ⁽¹⁰⁾		
Bovine Herpes Virus (BHV-1 and BHV-5)	Not detected	Not detected
Pestivirus (e.g. BVD-V)	Not detected	Not detected
Parainfluenza Virus Type 3 (PI-3)	Not detected	Not detected

Biological Performance			
Cell Growth MDCK (Canine Kidney Epithelial Cell Line)(11)	≥ 75	In testing	%
Cell Growth CHO (Chinese Hamster Ovary Cell Line)(11)	≥ 75	In testing	%

References/Methods

Freeze Point Depression (7)

(2) LAL-Test Membrane Filtration (total plate count, aerobic and anaerobic bacteria, yeasts and fungi)

Photometry

(3) Cytopathic Effect in Cell Culture and Immunofluorescence (4) Hexokinase Assay (10)

Biuret Colorimetry Cell Promoting Assay (5)(11)Nephelometry (6)

Intended Use For research use only. June 25, 2021 Date

This document has been created digitally and is valid without a signature.