

Certificate of Analysis

Product Fetal Bovine Serum Dialyzed Collected in South America, sterile-filtered

Cat. No. AB-FBS-DIA-0500, AB-FBS-DIA-1050, AB-FBS-DIA-0050

Lot. No.AB20214165Expiry Date03/2026OriginSouth AmericaManufacturing Date03/2021Storage≤-15°C

Quality Profile	Specification	Result	Unit
Appearance	Clear amber	Pass	
рН	6.8 – 8.2	7.38	
Osmolality ⁽¹⁾	280 – 340	281	mOsm/kg
Endotoxin ⁽²⁾	As reported	1.25	EU/mI
Hemoglobin	≤ 25	14.99	mg/dl
Glucose	≤ 0.28	0.28	mM
Calcium	As reported	0.59	mM
Phosphorous	As reported	0.35	mM
GOT/AST	As reported	44	U/I
LDH	As reported	2500	U/I

Protein Profile			
Total Protein ⁽³⁾	3.0 – 4.5	3.60	g/dl
Albumin ⁽⁴⁾	As reported	1.75	g/dl
IgG	≤ 400.0	191.06	mg/l

Sterility Testing			
Bacteria and Fungi	Not detected	Not detected	
Mycoplasma ⁽⁵⁾	Not detected	Not detected	

Virus Testing ⁽⁶⁾		
Bovine Herpes Virus (BHV-1 and BHV-5)	Not detected	Not detected
Pestivirus (e.g. BVD-V)	Not detected	Not detected
Parainfluenza Virus Type 3 (PI-3)	Not detected	Not detected

Biological Performance			
Cell Growth L929 (Murine Fibroblast Cell Line)(7)	≥ 75.0	101.0	%
Cell Growth HELA (Human Epithelial Cell Line)(7)	≥ 75.0	98.0	%

References/Methods

(1) Freeze Point Depression (5) Cell Culture

(2) Chromokinetic Test (6) Cytopathic Effect in Cell Culture and ELISA

(3) Biuret Colorimetry (7) Growth Promoting Assay

(4) Immunoturbidimetry

Intended Use For research use only.

Date June 01, 2021

This document has been created digitally and is valid without a signature.