

Certificate of Analysis

Product	Fetal Bovine Serum Collected in South America, sterile-filtered			
Cat. No.	AB-FBS-0500 S(500 ml), AB-FBS-1050S (500 ml), AB-FBS-0050S (50 ml)			
Lot. No.	AB20214448	Expiry Date	11/2026	
Origin Storage	South America ≤-15°C	Manufacturing Date	11/2021	

Quality Profile	Specification	Result	Unit
Appearance	Clear amber	Pass	
рН	6.8 - 8.2	7.24	
Osmolality ⁽¹⁾	280 - 340	309	mOsm/kg
Endotoxin ⁽²⁾	As reported	2.86	EU/ml
Hemoglobin ⁽³⁾	≤ 25.0	9.30	mg/dl
Glucose ⁽⁴⁾	As reported	3.55	mM
Calcium ⁽³⁾	As reported	2.87	mM
Phosphorous ⁽³⁾	As reported	2.62	mM
GOT/AST ⁽³⁾	As reported	31	U/I
LDH ⁽³⁾	As reported	506	U/I

Protein Profile			
Total Protein ⁽⁵⁾	3.0 – 4.5	3.9	g/dl
Albumin ⁽⁶⁾	As reported	2.71	g/dl
IgG ⁽⁷⁾	≤ 250	92.0	mg/l

Sterility Testing			
Bacteria and Fungi ⁽⁸⁾	Not detected	Not detected	
Mycoplasma ⁽⁹⁾	Not detected	Not detected	

Virus Testing ⁽¹⁰⁾		
Bovine Herpes Virus (BHV-1 and BHV-5)	Not detected	Not detected
Pestivirus (e.g. BVD-V)	Not detected	Not detected
Parainfluenza Virus Type 3 (PI-3)	Not detected	Not detected

Biological Performance			
Cell Growth MDCK (Canine Kidney Epithelial Cell Line)(11)	≥ 75.0	94.0	%
Cell Growth CHO (Chinese Hamster Ovary Cell Line)(11)	≥ 75.0	103.0	%

Refer	ences/Meth	ods		
(1)			(7)	ELISA
(2)	LAL-Test	·	(8)	Membrane Filtration (total plate count, aerobic and anaerobic bacteria, yeasts and fungi)
(3)	(3) Photometry		(9)	PCR
(4)			(10)	Cytopathic Effect in Cell Culture and Immunofluorescence
(5)	(5) Biuret Colorimetry		(11)	Cell Promoting Assay
(6) Nephelometry				
Inton	ded Use	For research use only		
	ueu use	For research use only.		
Date		January 13, 2022		
This document has been created digitally and is valid without a signature.				