

Certificate of Analysis

Product Fetal Bovine Serum, Collected in South America, sterile-filtered

Cat. No. AB-FBS0500S, AB-FBS-0050S, AB-FBS-1050S

Lot. No. AB20214403 Expiry Date 09/2026
Origin South America Manufacturing Date 09/2021

Storage ≤-15°C

Quality Profile	Specification	Result	Unit
Appearance	Clear amber	Pass	
рН	6.8 – 8.2	7.35	
Osmolality ⁽¹⁾	280 – 340	312	mOsm/kg
Endotoxin ⁽²⁾	As reported	1.93	EU/mI
Hemoglobin ⁽³⁾	≤ 25	10.70	mg/dl
Glucose ⁽⁴⁾	As reported	3.66	mM
Calcium ⁽³⁾	As reported	2.94	mM
Phosphorous ⁽³⁾	As reported	2.62	mM
GOT/AST ⁽³⁾	As reported	30	U/I
LDH ⁽³⁾	As reported	501	U/I

Protein Profile			
Total Protein ⁽⁵⁾	3.0 – 4.5	3.6	g/dl
Albumin ⁽⁶⁾	As reported	2.1	g/dl
IgG ⁽⁷⁾	≤ 250	83.31	mg/l

Sterility Testing			
Bacteria and Fungi ⁽⁸⁾	Not detected	Not detected	
Mycoplasma ⁽⁹⁾	Not detected	Not detected	

Virus Testing ⁽¹⁰⁾		
Bovine Herpes Virus (BHV-1 and BHV-5)	Not detected	Not detected
Pestivirus (e.g. BVD-V)	Not detected	Not detected
Parainfluenza Virus Type 3 (PI-3)	Not detected	Not detected

Biological Performance			
Cell Growth MDCK (Canine Kidney Epithelial Cell Line)(11)	≥ 75	100.0	%
Cell Growth CHO (Chinese Hamster Ovary Cell Line)(11)	≥ 75	99.0	%

References/Methods

) Freeze Point Depression (7) ELI

(2) LAL-Test (8) Membrane Filtration (total plate count, aerobic and anaerobic bacteria, yeasts and fungi)

(3) Photometry (9) PCR

(4) Hexokinase Assay (10) Cytopathic Effect in Cell Culture and Immunofluorescence

(5) Biuret Colorimetry(6) Nephelometry(11) Cell Promoting Assay

Intended Use For research use only.

Date November 11, 2021

This document has been created digitally and is valid without a signature.