

Certificate of Analysis (Special)

Product Fetal Bovine Serum IgG , sterile-filtered

Cat. No. AB-FBS-IG-0500, AB-FBS-IG-1050

Lot. No. AB20214417 Storage \leq -15°C Manufacturing Date 2021.09 Expiry Date 2026.09

Quality Profile	Specification	Result	Unit
Appearance	Clear amber	Pass	
рН	6.8 – 8.2	7.46	
Osmolality ⁽¹⁾	280 – 340	320	mOsm/kg
Endotoxin ⁽²⁾	As reported	1.8	EU/ml
Hemoglobin	≤ 25	5.54	mg/dl
Glucose	As reported	2.12	mM
Calcium	As reported	3.07	mM
Phosphorous	As reported	2.53	mM
GOT/AST	As reported	38	U/I

Protein Profile			
Total Protein ⁽³⁾	3.0 – 4.5	3.3	g/dl
Albumin ⁽⁴⁾	As reported	2.0	g/dl
IgG	≤ 5.0	0.14	μg/ml

Sterility Testing			
Bacteria and Fungi	Not detected	Not detected	
Mycoplasma ⁽⁵⁾	Not detected	Not detected	

Virus Testing ⁽⁶⁾			
Bovine Herpes Virus (BHV-1 and BHV-5)	Not detected	Not detected	
Pestivirus (e.g. BVD-V)	Not detected	Not detected	
Parainfluenza Virus Type 3 (PI -3)	Not detected	Not detected	

Biological Performance			
Cell Growth 929 (Murine Fibroblast Cell Line) ⁽⁷⁾	≥ 75.0	100.0	%
Cell GrowthHELA (Human Epithelial Cell Line) (7)	≥ 75.0	100.0	%

References/Methods

(1) Freeze Point Depression

(5) Cell Culture

(2) Chromokinetic Test

(6) Cytopathic Effect in Cell Culture and ELISA

(3) Biuret Colorimetry

(7) Growth Promoting Assay

(4) Immunoturbidimetry

Intended Use For research use only.

This document has been created digitally and is valid without a signature.